



NO ANIMAL NO HUMAN JUST TECHNOLOGY

From Design to Delivery

The Fastest AI-powered Organ-on-chip solution



Aligned with FDA's 2025
Alternative Testing Roadmap



WHERE BIOLOGY MEETS ENGINEERING: THE FUTURE OF IN VITRO MODELING



DRIVEN BY Taiwan's Semiconductor BUILD ORGANS AT SPEED



Biomimetic Organ-on-Chip Platform

- Validated in **10+ patient-derived** organ models, including airway, liver, intestine, and cancer.
- Supports **10+ biological indicators**, such as barrier integrity, mucus production, inflammatory response, and drug-induced toxicity.

Dynamic Perfusion System for Physiological Simulation

- Maintains dynamic culture with programmable flow and physiological shear stress.
- Supports cell differentiation and sustained functionality **across the 1-28 day culture** period.




Aerosol Exposure Module for Inhalation

- Automates 8-chip culture with **air, liquid, and aerosol delivery under air-liquid conditions**.
- Enables inhalation studies on drug deposition, ciliary motion, mucus clearance, and particle behavior for respiratory therapy testing.




Technology Integration Platform


The Fastest Path From Organ Modeling To Toxicology Decisions



Multi-organ type validation



Powering the Fastest Toxicity Profiling




By leveraging our extensive library of organ models, we enables **early toxicity detection** within just 6 hours— **faster** than traditional method.

Unlocking Clinical Insights

An AI agent built for organ-on-chip platforms, streamlining the path from design to clinical decisions

Concentration-Response Relationship

Figure 1: The relationship between different concentrations of Survanta and the measured %Area of SP-B.



Effect of Nebulized Survanta on SP-B Distribution in A549 Cells

The effect of nebulized Survanta on SP-B distribution in A549 cells was assessed by measuring the percentage area of SP-B immunofluorescence. As shown in Table 1, control cells exhibited a baseline SP-B distribution of $17.20 \pm 2.17\%$ area. Treatment with nebulized Survanta at 3.125 mg/ml resulted in a slight decrease in SP-B distribution to $15.66 \pm 2.54\%$ area, although this change was not statistically significant compared to control ($p > 0.05$).

Table 1: Percentage Area of SP-B Fluorescence in A549 Cells Treated with Different Concentrations of Nebulized Survanta



Doctor

Preclinical to Clinical Transition

Safety

Toxicity

Speciality: preclinical simulation and AI based clinical translation

Our Global Ecosystem

Accelerating Every Stage of Drug Development

NAR Labs
國家實驗研究院
National Applied Research Laboratories

TSRI **NAR Labs** 國家實驗研究院
台灣半導體研究中心
Taiwan Semiconductor Research Institute

DARWIN

LARGAN

臺中榮民總醫院
Taichung Veterans General Hospital



秀傳醫療體系
Show Chwan Healthcare System

馬偕紀念醫院
Mackay Memorial Hospital

國防醫學院
NATIONAL DEFENSE MEDICAL CENTER

台大醫院
NTUH

長庚紀念醫院
Chang Gung Memorial Hospital

天主教永和耕莘醫院
Yonghe Cardinal Tien Hospital

Cellentia

tlc
台灣微脂體

MOLECULAR DEVICES

PHISON

比翼生醫創投
BE HEALTH VENTURES



NAIP 北美智權
North America Intellectual Property Co.

CONTACT US

Phone: +886 3 6126 258

Web: www.anivance.io

Email: info@anivance.io

Anivance AI's platform supports the transition from discovery to regulatory submission, enabling human-relevant preclinical evaluation aligned with IND/CMC and FDA's New Approach Methodologies (NAMs).

Regional Offices

Taiwan / Headquarters: Rm. 6, 2F., No. 150, Sec. 2, Shengyi Rd., Zhubei City, Hsinchu County 302058, Taiwan.

Check our website for a current listing of worldwide distributors.